

Power: MAX 5W

Voltage: 3V

Battery: 450 mAh

Frequency: 8000 ($\pm 10\%$) rpm

Net weight: 100 g

Size: 140 x 183 x 66 mm

Standby time: 45 days

ORION



The Spectra Sculpt® kit contains: ORION Device. | USB Cable | Microfiber Cloth. | 3 pairs - 6 Multiplex Conductive Pads. | User Manual

Visit us at SPECTRASCULPT.com
to learn more about Quantum EMS Technology.

Registered trademark, Patented
CE, FCC & ROHS & FDA Certified

Executive Standards: GBT36419-2018

SPECTRA SCULPT, Brandstrasse 31, 8952 Schlieren, Switzerland



ENGLISH

Spectra Sculpt® had a vision for effective skincare and the cellular aging process while nurturing and providing deep cellular rejuvenation, division, and repair with exceptional results.

The precise frequency and impulse from our intelligent built-in program, trigger the skin's own ability to regenerate itself, that further stimulates additional ones, eventually creating an effective loop of rejuvenation across skin layers and muscle tissue to fight the appearance of all signs of aging. Spectra Sculpt® ORION is a non-invasive at-home facial toning device intended for cosmetic use.

ORION comes with an intelligent built-in program that runs for 15 minutes, with perfectly calibrated Quantum EMS, and allows users to adjust the output level for personalized comfort.

-
- + Reducing Nasolabial Folds
 - + Lifting Cheeks
 - + Defining Jawline
 - + Tightening Jowls
 - + Lifting Eyebrows & Eyelids
 - + Fuller & Plumper Lips
 - + Uplifted & Raised Corner of the Mouth
 - + Reduction of Overall Puffiness
 - + Tightening and Toning the Periorbital Eye Area
 - + Reducing Fine Lines and Wrinkles noticeably
 - + Diminish Cellulite
 - + Muscle recovery

WARRANTY

WHAT WE'LL NEED:

Full Name

Shipping address

Phone number

Product to be returned

Serial number – located on the back of the box. Scratch security seal to reveal serial number. We kindly suggest to note your serial no. here for your record. SERIAL No.: _____

Description of the issue you are experiencing with your device.

Proof of purchase – please upload a copy or image of your proof of purchase from an authorized seller.

Email the above to
customercare@spectrasculpt.com for verification.

WARRANTY

SPECTRA SCULPT® manufactures its products to the highest standards. Each device undergoes rigorous design, safety, inspection, and quality tests for every component before they reach you at the consumer level.

Every SPECTRA SCULPT® device comes with a serial number to certificate authenticity. We kindly advise you to register your device with your serial no. at www.spectrasculpt.com/register to ensure 100 % coverage.

If your SPECTRA SCULPT® device does not operate correctly or got damaged during shipping, please follow our warranty & return policy below:

SPECTRA SCULPT® will replace its devices qualified for warranty. Every purchased device includes one (1) year of warranty for any manufacturer-related problems from normal use of the device and is not refreshed with replacement devices. Warranties are only valid if returned within the warranty period, have a serial number and a receipt. The warranty will be considered void if the product has been subjected to an accident, misuse, improper repair, or any other use not listed in the printed instructions.

If you did not purchase your device through our website, please contact the dealer directly.

PRODUCT INSTRUCTIONS

Please be sure to fully charge ORION before use. This may take up to 8 H. The standby time after a full charge is 45 days.

Connect one end of the USB cable to ORION and the other end to a USB charger port. The red light on ORION will flash while charging.

Always clean the device with the dry microfiber cloth before and after use.

Use ORION for 7-9 minutes - max. 15 minutes - 4-5 times a week for a minimum of 8 weeks. After this period, use 1-3 times a week to maintain the results or as required. Mature users, continue to see improvements for up to 4 months.

Reduction of cellulite and tightening of loose skin around the knee and legs can be seen in as little as 3 weeks. Continue to use for max. 15 minutes for a minimum of 6 weeks for continuous improvement.

Use after a long day and/or heavy legs - wait for a minimum of 1 H after a workout before use for recovery.

Use before any event for an immediate lifting effector or quickly remove morning puffiness or built-up fluid for 3-5 minutes.

To familiarize yourself with the quantum EMS sensation for the first time, place ORION around your head and hold one massage head a few millimeters away. Be sure not to touch the metal part with your fingers. Turn ORION on and slowly place the second massage head on your face.

FACE & LEG PLACEMENT



The conductors can be placed anywhere in this area and do **not** need to be moved outside this region to target all facial areas including the neck.

General cellular repair | Level: 1 & 2

Circulation / intense treatment | Level: 3

LEVEL 1

Low

LEVEL 2

Moderate

LEVEL 3

Max. + 8000 rpm

Increase treatment time and level slowly and stop where it feels comfortable for you. Your skin might only tolerate a certain amount of electrical impulses depending on the density of your tissue.



The conductors can be placed anywhere in this area - aligned to the front or the back. Do not place it directly on the knee or the groin.

Cellulite | Level: 3

Regeneration / recovery | Level 1 & 2

GENERAL

Your Spectra Sculpt® device is designed for cosmetic use only and for individuals in good health. Please read the instruction manual carefully before using this device.

ADVERSE REACTIONS

Some users with sensitive skin may experience skin irritation or hypersensitivity due to electrical stimulation. The irritation can usually be avoided by lowering the intensity level and treatment time. Allergic reactions, skin irritation, or burns beneath the electrodes have been reported with powered muscle stimulators. If you experience any allergic reaction, irritation or burning, discontinue use.

VERY IMPORTANT

Any harmful consequence resulting from misuse, disregarding safety instruction, connection to improper voltage sources, dirty conductors, or any other improper applications is neither the responsibility of Spectra Sculpt nor its affiliates or distributor.

WARNING

Handle ORION with care. It contains sensitive electronic components that may get damaged.

Do not drop ORION.

Do not use in environment of excessive humidity or dense dust

Do not use in environments above 40° or below 5° degrees Celsius.

Do not charge if the USB cable is damaged.

Keep away from water.

Do not use ORION if it has been in contact with water.

Do not clean ORION with water.

Keep away from heat.

Keep away from metal.

Keep away from children.

Keep away from flammable substances.

Do not disassemble ORION.

Do not use for any purpose other than those specified in the manual.

Do not use it while operating any machinery.

Do not place magnets next to ORION.

Caution when there is skin that lacks normal sensation.



Gently wash your face. Do not use any oil-based products or heavy skincare before use as they will act as a barrier.

Users with facial hair or excessive amount of peach fuzz should shave before treatment, as hair can interfere with the conductivity and can cause irritations.

Remove any jewelry and/or piercings before use. Use the multiplex conductive gel pads underneath the conductors. Pads are 10-15 times reusable.



Clean the device with the dry microfiber cloth. Place the multiplex conductive gel pad on the conductors.

Place ORION around your head and make sure both massage heads are placed well on your cheeks.

Adjust the handles to your comfort

Alternative, ORION can be placed on the lower part of your thigh, to help fight cellulite and sagging skin.



Enjoy your treatment 7-9 minutes - maximum 15 minutes. Modify time to your comfort. Tolerance may vary from user to user. Treatment should never feel uncomfortable!

IMPORTANT - Start slow. If you are new to using ORION, begin with 3-5 min. on level 1 and work your way up to avoid overstimulation and/or irritation.

MODE OF OPERATION

To turn it ON, long press the button for 2 sec. You will hear a sound signal as confirmation.

To switch to level 2, quickly press the button for 1 sec. You will hear 2 sound signals as confirmation.

To switch from level 2 to 3, quickly press again for 1 sec and you will hear 3 sound signals as confirmation.

To turn it OFF, press the button for 3 sec. The red light will be turned off, indicating a successful power off.

When fully charged, the red light will be permanently on and stop flashing. In case of low battery, ORION will send out an alarm signal.

After 15 minutes, ORION will automatically turn OFF.

Do not exceed 15 minutes in 24 H, your body needs time to build its new structure.



SAFETY PRECAUTIONS

DO NOT USE if the following conditions apply: Patient with heart disease, cardiac pacemaker, any electronic implanted devices, high fever, cancer, a recent history of cancer, tumor, infectious disease, recent skin wound, any kind of skin disease (acne, subcutaneous inflammations), acute illness, pregnancy, sunburn, vasodilation, metal implants, epilepsy, and simultaneous use of medical electrical devices or any other electrical device (e.g., ear-pods).

Do not use if you have acne or are prone to inflammation on the skin or any skin conditions.

Do not use without conductive gel pads if you have sensitive skin.

Do not use if you have any dental metal implants or metal fillings.

Do not place conductors on your neck or move them up to your eyes / forehead.

Do not use if you have a history of epilepsy or seizures.

Do not use if you recently had plastic surgery or other surgeries.

Do not use if the skin is freshly injured or sunburned.

Stop using immediately if you get persistent headaches, feel irritated, dizzy, pain, or hot. Treatment should never be unpleasant.

Do not touch the metal part with your fingers while in use.

The long-term effects of stimulation are not known.